

Oticon Hearing Instruments and DAI/FM compatibility Overview 2006

910 53 310 0605.06 printed in Denmark

Hearing Aid	DAI Adaptor		FM Adaptor	FM receiver
	Low Level Input	High Level Input		
	Bodyworn FM systems e.g. Solaris	TV, stereo, MP3 player, CD player, Computer	Adaptor for AMIGO receivers and other ear-level plug-in receivers	Universal and dedicated AMIGO receivers
Oticon Syncro, Oticon Safran, Oticon Tego Pro, Oticon Tego 	AP 800 399-51-480-00* Transparent  * includes tool for adaptor lock * input levels are set manually		FM 8 399-51-470-00* Transparent  * includes tool for adaptor lock	R1/R2 Silver 
Adapto, Gaia, Atlas Plus Direct, GO Direct DigiFocus II (all V2*)  <small>* if not marked V2 on battery door, aid must be sent to Oticon for modification.</small>	AP 501 399-51-440-03* Transparent  * includes tool for making hole in battery door	AP 502 399-51-441-04* Transparent  * includes tool for making hole in battery door	FM 5 399-51-170-07* Transparent  * includes tool for making hole in battery door	R1/R2 Silver 
DigiFocus II Super Power, 300P 	AP 201 399-50-360-05 Beige 	AP 202 399-50-370-07 Beige 	FM 3 399-51-180-09* Transparent  * includes tool for making hole in battery door	R1/R2 Silver 

Hearing Aid	DAI Adaptor		FM Adaptor	FM receiver
	Low Level Input	High Level Input		
	Bodyworn FM systems e.g. Solaris	TV, stereo, MP3 player, CD player, Computer	Adaptor for AMIGO receivers and other ear-level plug-in receivers	Universal and dedicated AMIGO receivers
Personic 	AP 401 399-50-276-02 Transparent 	AP 402 399-50-286-04 Transparent 	FM 4 399-51-130-08* Transparent  * includes tool for making hole in battery door	R1/R2 Silver 
Atlas, Atlas Plus, GO, (all except Direct) Ergo, Swift 	AP 600 399-50-460-06* Transparent  * includes tool for adaptor lock * input levels are set manually		FM 6 399-50-450-04* Transparent  * includes tool for adaptor lock	R1/R2 Silver 
Sumo DM, Sumo XP, Sumo E, Swift 120+ 	AP 700 399-50-900-08* Grey  * input levels are set manually		FM 7 399-50-920-02 Grey 	R1/R2/R7 Silver 
380P, 390PL 	AP 301 399-50-340-00 Beige 	AP 302 399-50-350-03 Beige 	FM 3 399-51-180-09* Transparent  * includes tool for making hole in battery door	R1/R2 Silver 
CI Processors from various manufacturers 			FM-CI Cochlear implant adaptor 140-08-110-00 	R1/R2 Silver 



FM receiver order information:	R1:	169-187 MHz	216-217 MHz
R1:	140-06-000-00	140-06-100-00	
R2:	140-06-200-00	140-06-300-00	
R7:	140-06-600-00	140-06-700-00	